


# Bottles for Injectables - PP

50ml, 100ml, 250ml

**PLASTIC**  
Compatible  
with the  
environment



Ref. n°	1839	1840	1697C
Volume (ml)	50	100	250
Material	PP	PP	PP
Colour	Natural	Natural	Natural
Dimension A (mm)	64.6	87.6	121.5
Dimension B (mm)	41	61	77.5
Dimension C (mm)	Ø43	Ø49	Ø66
Neck opening (mm)	Ø12.7	Ø12.7	Ø12.7
Weight (gr)	10.1	14.6	30
Packing (pcs/brick)	55	41	35
Packing (pcs/carton)	605	492	210
Volume of carton (m3)	0.094	0.144	0.144
Carton dimension l x w x h (mm)	600x400x391	800x400x450	800x400x450
Carton weight filled (kg)	6	7	73
Cartons/EUR pallet	20	15	15
Quantity/EUR pallet	12,100	7,380	3,150
Sterilization	Autoclavable		
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6		

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.