


Vials for Injectables Polypropylene (PP)

50ml, 100ml



| Ref. n° | 1917 | 1840 |
|------------------------------|---|---------------|
| Volume (ml) | 50 | 100 |
| Material | PP | PP |
| Colour | Natural/Amber | Natural/Amber |
| Dimension A (mm) | 71.5 | 87.6 |
| Dimension B (mm) | 41 | 61 |
| Dimension C (mm) | Ø43 | Ø49 |
| Neck opening (mm) | Ø12.7 | Ø12.7 |
| Weight (gr) | 11.0 | 14.6 |
| Brimful Volume (ml) | 59.0 | 111.0 |
| Packing (pcs/brick) | 55 | 41 |
| Packing (pcs/carton) | 440 | 492 |
| Volume of carton (m3) | 0.094 | 0.144 |
| Carton dimension LxWxH (mm) | 600x400x391 | 800x400x450 |
| Carton weight filled (kg) | 6 | 7 |
| Cartons/EUR pallet | 20 | 15 |
| Quantity/EUR pallet | 8,800 | 7,380 |
| Sterilization | Autoclavable | |
| Rubber disc or stopper - Cap | The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6 | |

2021-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.