

Clear Vial™

5 ml, 10 ml, 20 ml, 50 ml, 100 ml

PLASTIC
Compatible
with the
environment



Ref. no.	1875	1876	1877	1879	1880
Volume	5 ml	10 ml	20 ml	50 ml	100 ml
Material	COC/COP	COC/COP	COC/COP	COC/COP	COC/COP
Colour	Natural/ Amber	Natural/ Amber	Natural/ Amber	Natural/ Amber	Natural/ Amber
Dimension A (mm)	40	45	55	73	94.5
Dimension B (mm)	Ø22	Ø24	Ø30	Ø42.5	Ø52.6
Dimension C (mm)	Ø20	Ø20	Ø20	Ø20	Ø20
Weight (gr)	2.74	3.92	5.48	10.1	14.8
Brimful Volume	9.0	14.0	25.0	72.0	139.0
Packing (pcs/brick)	230	162	135	55	36
Packing (pcs/carton)	4,140	2,592	1,620	770	432
Volume of carton (m3)	0.094	0.094	0.094	0.094	0.144
Carton dimension LxWxH (mm)	600x400x391	600x400x391	600x400x391	800x400x450	800x400x450
Carton weight filled (kg)	12.3	11.2	8.5	6.0	7.0
Cartons/EUR pallet	20	20	20	15	15
Quantity/EUR pallet	82,800	51,840	32,400	11,550	6,480
Sterilization	Irradiation/Autoclaving				
Rubber disc or stopper and cap - The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6					

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.