


Bottles for Injectables PE

250ml, 500ml, 1000ml

PLASTIC
Compatible
with the
environment



Ref. no	1697	1698	1699
Volume (ml)	250	500	1000
Material	PE-HD	PE-HD	PE-HD
Colour	Natural	Natural	Natural
Dimension A (mm)	121.8	157.6	203.3
Dimension B (mm)	77.5	106.2	141.8
Dimension C (mm)	Ø66	Ø78	Ø95
Neck opening (mm)	Ø22.5	Ø22.5	Ø22.5
Weight (gr)	28	38	60
Brimful Volume	274.0	535.0	1,043.0
Packing (pcs/brick)	35	27	16
Packing (pcs/carton)	210	135	64
Volume of carton (m3)	0.144	0.144	0.144
Carton dimension l x w x h (mm)	800x400x450	800x400x450	800x400x450
Carton weight filled (kg)	6.5	6.6	5.3
Axial strength (filled 20°C) - (kg)	100	100	100
Cartons/EUR pallet	15	15	15
Quantity/EUR pallet	3,150	2,025	960
Sterilization	Beta ray		
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8536-2 and aluminium caps according to ISO 8536-3/7		

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.