


Vials for Injectables PP

3ml, 5ml, 10ml, 20ml, 25ml

PLASTIC
Compatible
with the
environment



Ref. n°	1834	1835	1836	1837	1838
Volume (ml)	3	5	10	20	25
Material	PP	PP	PP	PP	PP
Colour	Natural / Amber	Natural / Amber	Natural / Amber	Natural / Amber	Natural / Amber
Dimension A (mm)	25.9	37.9	45.6	48.6	53.5
Dimension B (mm)	14	26	32	29	38
Dimension C (mm)	Ø22	Ø22	Ø26	Ø32	Ø33
Neck opening (mm)	Ø12.6	Ø12.6	Ø12.6	Ø12.6	Ø12.6
Weight (gr)	2.05	2.80	3.92	5.50	6.25
Brimful Volume	5.5	9.0	15.3	23.5	29.6
Packing (pcs/brick)	230	230	162	98	85
Packing (pcs/carton)	5,760	4,140	2,592	1,372	1,105
Volume of carton (m3)	0.094	0.094	0.094	0.094	0.094
Carton dimension lwxhxh (mm)	600x400x391	600x400x391	600x400x391	600x400x391	600x400x391
Carton weight filled (kg)	12.5	12.3	11.2	8.5	7.9
Cartons/EUR pallet	20	20	20	20	20
Quantity/EUR pallet	115,000	82,800	51,840	27,440	22,100
Sterilization	Autoclavable				
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-316				

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.