


Bottles for Injectables - PP

50ml, 100ml

PLASTIC
Compatible
with the
environment



Ref. n°	1917	1840
Volume (ml)	50	100
Material	PP	PP
Colour	Natural/Amber	Natural/Amber
Dimension A (mm)	71.5	87.6
Dimension B (mm)	41	61
Dimension C (mm)	Ø43	Ø49
Neck opening (mm)	Ø12.7	Ø12.7
Weight (gr)	11.1	14.6
Brimful Volume (ml)	58.0	111.0
Packing (pcs/brick)	-	41
Packing (pcs/carton)	440	492
Volume of carton (m3)	0.094	0.144
Carton dimension l x w x h (mm)	600x400x391	800x400x450
Carton weight filled (kg)	6	7
Cartons/EUR pallet	20	15
Quantity/EUR pallet	8,800	7,380
Sterilization	Autoclavable	
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6	

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.