

# Clear Vial™

5 ml, 10 ml, 20 ml, 50 ml, 100 ml



Ref. no.	1875	1876	1877	1879	1880
<b>Volume</b>	5 ml	10 ml	20 ml	50 ml	100 ml
<b>Material</b>	COC	COC	COC	COC	COC
<b>Colour</b>	Natural	Natural	Natural	Natural	Natural
<b>Dimension A (mm)</b>	40	45	55	73	94.5
<b>Dimension B (mm)</b>	Ø22	Ø24	Ø30	Ø42.5	Ø51.6
<b>Dimension C (mm)</b>	Ø20	Ø20	Ø20	Ø20	Ø20
<b>Weight (gr)</b>	2.9	3.5	5.1	12.5	19.5
<b>Brimful Volume (ml)</b>	10	14	27	70	134
<b>Packing (pcs/brick)</b>	230	162	135	55	36
<b>Packing (pcs/carton)</b>	4,140	2,592	1,620	770	432
<b>Volume of carton (m3)</b>	0.094	0.094	0.094	0.144	0.144
<b>Carton dimension LxWxH (mm)</b>	600x400x391	600x400x391	600x400x391	800x400x450	800x400x450
<b>Carton weight filled (kg)</b>	12.8	9.9	9.0	10.4	9.2
<b>Cartons/EUR pallet</b>	20	20	20	15	15
<b>Quantity/EUR pallet</b>	82,800	51,840	32,400	11,550	6,480
<b>Sterilization</b>	Irradiation/Autoclaving				
<b>Rubber disc or stopper - Cap</b>	The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6				

2021-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.