


# Vials for Injectables Polyethylene (HDPE)

50 mL, 100 mL



Ref. no.	1719	1680
Volume (mL)	50	100
Material	HDPE	HDPE
Colour	Natural	Natural
Dimension A (mm)	66.5	82.1
Dimension B (mm)	35.0	45.0
Dimension C (mm)	Ø43	Ø52
Neck opening (mm)	Ø12.7	Ø12.7
Weight (gr)	11.0	14.0
Brimful Volume (mL)	55	106
Packing (pcs/pack unit)	121	96
Packing (pcs/carton)	605	288
Volume of carton (m <sup>3</sup> )	0.094	0.094
Carton dimension LxWxH (mm)	600x400x391	600x400x391
Carton weight filled (kg)	7.6	5.0
Cartons/EUR pallet	20	20
Quantity/EUR pallet	12,100	5,760
Sterilization	Irradiation	
Rubber disc or stopper - Cap	The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6	

2024-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.