


# Bottles for Injectables Polypropylene (PP)

100ml, 250ml, 500ml, 1000ml



Ref. n°	1894	1997	1698B	1699B
Volume (ml)	100	250	500	1000
Material	PP	PP	PP	PP
Colour	Natural/Amber	Natural/Amber	Natural/Amber	Natural/Amber
Dimension A (mm)	98.5	121.8	157.6	203.3
Dimension B (mm)	64.4	77.5	106.2	141.8
Dimension C (mm)	Ø49	Ø66	Ø78	Ø95
Neck opening (mm)	Ø22.5	Ø22.5	Ø22.5	Ø22.5
Weight (gr)	17	30	43	62
Brimful Volume (ml)	127.0	279.0	531.0	1053.0
Packing (pcs/brick)	36	35	27	16
Packing (pcs/carton)	432	210	135	64
Volume of carton (m3)	0.144	0.144	0.144	0.144
Carton dimension LxWxH (mm)	800x400x450	800x400x450	800x400x450	800x400x450
Carton weight filled (kg)	8.3	7.4	7.3	5.5
Cartons/EUR pallet	15	15	15	15
Quantity/EUR pallet	6,480	3,150	2,025	960
Sterilization	Autoclavable			
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8536-2 and aluminium caps according to ISO 8536-3/7			

2021-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.